

## Olmesartan Medoxomil Suspension 2 mg/mL

### INGREDIENTS:

Olmesartan medoxomil 20 mg tablet	12 tablets
Purified water USP	30 mL
Ora-Plus	30 mL
Ora-Sweet	QSAD: 120 mL

### EQUIPMENT AND SUPPLIES:

Graduated cylinder

### PREPARATION DETAILS:

1. Place tablets in an appropriately sized amber bottle.
2. Add purified water USP to amber bottle and allow to stand for a minimum of 5 minutes.
3. Shake the container for at least 1 minute and allow the suspension to stand for at least 1 minute.
4. Repeat shaking/standing process four additional times.
5. Add Ora-Plus and Ora-Sweet to the container and shake for at least 1 minute.

**Quality-Control Procedures** — Visually inspect for physical appearance of formulation and container closure integrity (no leakage, cracks in container, or improper seals).

**Labeling Requirements** — Extemporaneously compounded preparation. For oral use only. Store in refrigerator. Shake well before use.

**Storage Conditions/Stability** — Refrigerate. Stable for 28 days.

### STABILITY STUDY DETAILS:

**Study Container Type** — Amber polyethylene terephthalate (PET) bottle

**Referenced Manufacturers** — Olmesartan medoxomil tablets (Benicar, Daiichi Sankyo, Inc); purified water USP (not specified); Ora-Plus, Ora-Sweet (Paddock Laboratories, LLC).

**Stability-Indicating Study** — No

### REFERENCE

1. Benicar [prescribing information]. Basking Ridge, NJ: Daiichi Sankyo, Inc; 2019.