

Valsartan Suspension 4 mg/mL—Formulation 1

INGREDIENTS:

Valsartan 80 mg tablet	8 tablets
Ora-Plus	80 mL
Ora-Sweet SF	QSAD: 160 mL

EQUIPMENT AND SUPPLIES:

Graduated cylinder

PREPARATION DETAILS:

1. Place tablets in an appropriately sized amber bottle.
2. Add Ora-Plus to amber bottle.
3. Shake for at least 2 minutes, allow the suspension to stand for a minimum of 1 hour, and then shake for a minimum of 1 additional minute.
4. Add vehicle and shake well for at least 10 seconds to disperse contents.

Quality-Control Procedures — Visually inspect for physical appearance of formulation and container closure integrity (no leakage, cracks in container, or improper seals).

Labeling Requirements — Extemporaneously compounded preparation. For oral use only. Store at room temperature or refrigerate. Shake well before use.

Storage Conditions/Stability — Store at room temperature or refrigerate. Stable for 75 days when refrigerated and 30 days when stored at room temperature.

STABILITY STUDY DETAILS:

Study Container Type — Amber glass bottle with a child-resistant screw-cap closure

Referenced Manufacturers — Valsartan tablets (Diovan, Novartis Pharmaceuticals); Ora-Plus, Ora-Sweet SF (Paddock Laboratories, LLC).

Stability-Indicating Study — No

REFERENCE

1. Diovan [product information]. East Hanover, NJ: Novartis Pharmaceuticals; 2019.